

Use this guide to help select Spex products compatible to fit to the **Levo C3** powered standing wheelchair.

Step 1

Select a Compatible Back Support Size

Back Support	Width	Compatible Height*	
		Excluding Head Support Hardware	Including Stylo Slim Mounting Block
Spex (Classic, Manta, Vigour, SuperShape)	14 - 21"	14 - 25"	16 - 25"
Zygo Mid	14 - 20"	14 - 18"	16 - 18"


*When seat cushion height is the same as the Spex Classic Cushion.

Not compatible with Zygo Lo, Adapta², Adapta³, XLella, XLella² or Height-Adjustable Back Supports.

ⓘ Take Note of Compatibility Guide Assumptions

Spex and Manta Back Supports

The bottom edge of the back support is positioned flush with the seat cushion with the wheelchair in standing position.



Zygo Mid Back Supports

The bottom edge of the back support is positioned at least 3" above the cushion when the wheelchair is in seated position.

