

Anterior Trunk Positioning

Unique harness designs suit a wide scope of requirements for postural restraint and support.



FUNCTIONAL

- Pads follow natural contours of the body, with sleek design to minimise sweating
- Soft edge binding to reduce chafing, while heavy duty webbing will resist the twist
- Maintain and encourage spinal alignment over a stable pelvis

	Size	Code	Spec (1)	Spec (2)	Spec (3)	W/C Width	Image
H-HARNESS	XSML	1409-6640-017		27 cm		8" - 12"	
	SML	1409-6641-017		30 cm		11" - 15"	
	MED	1409-6642-017		36 cm		14" - 17"	
	LGE	1409-6643-017		42 cm		17" - 20"	
	XLGE	1409-6644-017		48 cm		20" +	
RETRACTOR HARNESS	SML	1409-6621-017		32 cm		11" - 15"	
	MED	1409-6622-017		38 cm		14" - 17"	
	LGE	1409-6623-017		44 cm		17" - 20"	
	XLGE	1409-6624-017		50 cm		20" +	
VEST HARNESS	XSML	1409-6610-017	15 cm	24 cm	24 cm	8" - 12"	
	SML	1409-6611-017	20 cm	32 cm	32 cm	11" - 15"	
	MED	1409-6612-017	26 cm	40 cm	40 cm	14" - 17"	
	LGE	1409-6613-017	30 cm	48 cm	48 cm	17" - 20"	
CENTRE-POINT HARNESS	XSML	1409-6660-017	23 cm	26 cm	22 cm	8" - 12"	
	SML	1409-6661-017	29 cm	36 cm	32 cm	11" - 15"	
	MED	1409-6662-017	34 cm	40 cm	36 cm	14" - 17"	
	LGE	1409-6663-017	40 cm	45 cm	42 cm	17" - 20"	
	XLGE	1409-6664-017	45 cm	50 cm	50 cm	20" +	

SPEX HARNESS GUIDES

Recommended for any harness choice to prevent shoulder straps falling off shoulders.

Spex Harness Guides: 1400-6692-000

Spex Compact Harness Guides: 1400-6693-000



Spex Harness Guides	Spex Compact Harness Guides
XLella	Adapta2
Zygo Mid	Adapta3
>14" width:	<13" width:
Vigour Hi	Vigour Hi
Manta	Manta
Classic	Classic
SuperShape	SuperShape

HARNESS OPTIONS



H-HARNESS

Minimalist option with good adjustment



RETRACTOR HARNESS

Encourages retraction of the shoulders



VEST HARNESS

- Greater coverage across the chest area where support is needed
- Low neck design for safety
- Stomach clearance comfort cutout



CENTREPOINT HARNESS

- Follows natural contours of the body and has clearance for items such as stomach feeding tubes
- Can be independently operated by the user